



MR17 NEUROSURGICAL TABLE

The new standard in intraoperative surgical tables

An integration of the Hill-Rom TruSystem® 7500 column,
and the MR Neuro OR Tabletop.



MR Neuro OR Tabletop



Hill-Rom TruSystem® 7500 column


FEATURES

- » MR conditional for safe intraoperative imaging
- » Multifunctional design with interchangeable tabletops is suitable for a wide range of applications and disciplines
- » Trendelenberg, reverse Trendelenberg, tilt, and longitudinal ranges for flexible positioning during surgery and scanning; extender board also available
- » Single break point in tabletop at the anatomical hip location for optimal patient positioning
- » Communicates with theatre controls to ensure correct configuration and orientation for safe intraoperative MRI

BENEFITS

- » Patients can be optimally positioned for surgical site access and intraoperative imaging
- » Designed to provide highly ergonomic working methods and to improve efficiency of surgical workflow
- » Precise, electric control for reliable and intuitive operation
- » Intuitive coupling point system enables tabletops to be exchanged with ease
- » Fully integrates with the IMRIS head fixation portfolio and other products, such as the MR Neuro Accessory Adapter

Technical Specifications for MR17 Neurosurgical Table

Description	 MR conditional (safe for use with 1.5T and 3T MR scanners)
Positioning Ranges	Trendelenburg: $\pm 45^\circ$ Tilt: $\pm 30^\circ$ Longitudinal slide: 472 mm Transverse slide: 250 mm (125 mm to the right and left from level position) Rotation: 350° Break adjustment: 0° to 30°
Height Range	748 mm to 1,298 mm 29.4" x 51.1"
Tabletop Dimensions	Length with tabletop fully extended = 3103 mm Length without extended tabletop (tabletop in level position) = 2631 mm Width of the back section without side rails = 450 mm Total height without pad = 383 mm
Maximum Patient Weight	180 kg 397 lbs
Power	External supply = 100V-230V ~ 50-60Hz / 240V ~ 50Hz Consumption = max. 700 VA

Ordering Information

Please email the purchase order directly to orders@imris.com.

Part number with operator manual:

119489-000 MR17 Neurosurgical Table - US and Canada

119490-000 MR17 Neurosurgical Table - Europe

Additional Products

120025-000 MR Neuro Accessory Adapter
(allows for use of non-MR conditional products)

For more information, contact IMRIS at +1 763-203-6300, or visit our website for additional product information, www.imris.com.



Table may be tilted during scanning



Patient's head below level of table