

HEAD FIXATION ACCESSORIES

MR conditional attachments for expanded versatility

For versatile head positioning and secure attachment of accessory systems.

LINKAGE EXTENSION

FEATURES

Optional linkage extension length +7.5cm (transitional member) available for flexible positioning of the head based on patient's size and anatomy.

BENEFITS

Allows for head positioning with translational movement in x, y, z directions, as well as rotational swivel, tilt, and pivot adjustability.

Linkage extension
(+7.5cm)



Rocker Arm
Accessory



HFD100 WITH ROCKER ARM
ACCESSORY FOR 4-PIN FIXATION

ROCKER ARM ACCESSORY


FEATURES

Converts from 3-pin to 4-pin head fixation; option to use inner or outer pin holes on each rocker arm, depending on skull anatomy.

BENEFITS

Distributes skull clamp pressure across two points on the side where the rocker arm is added.

Starburst
Adapter



STARBURST ADAPTER CONNECTED TO HFD100
WITH RETRACTOR SYSTEM ATTACHED

STARBURST ADAPTER

FEATURES




Specifically designed for attaching accessory systems to the head fixation device.

BENEFITS

Two industry-standard connections allow for secure attachment of other surgical accessory systems.

HFD100 HEAD FIXATION ACCESSORIES

TECHNICAL SPECIFICATIONS

	Linkage Extension	Rocker Arm Accessory	Starburst Adapter
General	For additional range when positioning head	Converts 3-pin to 4-pin head fixation	For attaching accessory systems
Sterility	Non-sterile	Non-sterile	Non-sterile
Safety	 MR Conditional	 MR Conditional	 MR Conditional

ORDERING INFORMATION

Please email the purchase order directly to orders@imris.com.

HFD100 ACCESSORIES:

Part number	Name
118492-000	HFD100 Linkage Extension, +7.5cm
117948-000	HFD100 Rocker Arm Accessory
119012-000	HFD100 Starburst Adapter Kit
802502-000	Navigation Adapter, Stryker
802503-000	Navigation Adapter, Medtronic
802504-000	Navigation Adapter, Brainlab