

The IMRIS InVision™ 3T Recharge Operating Suite

featuring the Siemens 3T Skyra Fit upgrade from Verio



 **IMRIS**

Confidence visualized

THE GOLD STANDARD IN NEUROSURGICAL INTRAOPERATIVE MRI

A POWERFUL UPGRADE TO AN IMRIS OPERATING SUITE

Upgrading an IMRIS Operating Suite with a Siemens Verio MRI to one that features a Siemens Skyra Fit with Biomatrix Technology can extend the useful life of an IMRIS system and enable IMRIS intraoperative MRI value to be maintained at a clinical center. With the InVision™ 3T Recharge Operating Suite, the clinical space is optimized to accommodate the Skyra Fit upgrade, and all major IMRIS system components. The InVision 3T Recharge brings the latest state-of-the-art technology to your current IMRIS Operating Suite, extending the life of your initial investment.

The upgrade provides the latest MRI technology, further improving diagnostic imaging quality. Workflows are streamlined through AI powered imaging protocols and automatic image reconstruction post processing as background activities, enabling increases in productivity while reducing scan times.¹




Image courtesy of Siemens Healthineers

INVISION 3T RECHARGE OPERATING SUITE UPGRADES

- Magnet Bay Doors
- ORT-400 Operating Room Table
- Magnet Mover and Rail System
- Head Fixation Devices
- RF Shielding Modifications
- AP5 Room Control
- Collision Detection System
- On-Site Clinical Workflow & Safety Training Program

¹ <https://www.siemens-healthineers.com/en-us>

NEW FEATURES WITH THE IMRIS INVISION™ 3T RECHARGE MAGNET





THE INVISION OPERATING SUITE FROM IMRIS

MAGNETIC PARAMETERS

Operating field strength	3 Tesla
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogen)	6,400 kg \pm 100 kg
Magnet length	1,630 mm
Open Bore design	700 mm \pm 1%
System length cover to cover	1,730 mm

IMRIS

1230 Chaska Creek Way, Suite 100
Chaska, MN 55318

Main Line: +1 763-203-6300

Customer Support: +1 866-475-0525

Email: customersupport@imris.com

Sales: +1 763-203-6407

Email: sales@imris.com



Confidence visualized